

floorline® 9.5|80

A new concept for barrier-free care

Perfectly suited for care.

.bock®



New dimensions to adjustable height

– ideal for residents and care givers

The floorline® 9.5|80 can be lowered extremely close to the ground: in the ultra-low position of only 9.5 cm above the floor, it provides protection from falls while still giving the residents freedom to move. Moreover, the residents' ability to simply reach out and touch the floor helps them to feel safe and secure.

The floorline® 9.5|80 also makes care givers' everyday activities easier: in the ultra-low position, there is no need for otherwise common, additional safety measures such as laying out mats. What's more, this specially designed bed allows residents to be more mobile by facilitating the process of getting in and out of the bed without assistance.

// Reliable protection from falls

// Sliding instead of falling

// Sleepers can lie without restriction and how they wish

// Healthy and back friendly work



The barrier-free bed – for independent mobility and autonomy

The barrier-free lowered position enhances the individual feeling of safety and security for residents and allows for independent mobility –both when leaving the bed and when returning. On all fours, this is even possible without help for people with limited mobility.

Using the hand control, residents can select the comfortable bed height of 37.5 cm for climbing in and out of bed.

With the floorline[®] 9.5|80, measures to protect against falls that limit the residents' freedom of movement can generally be done without – in accordance with the Werdenfelser Approach. Moreover, nursing care beds from bock support the requirements of the future expert standards of mobility. If desired, various types of side rails offer support in helping residents out of bed, and additional safety, security and stability.

- // Convenient height for getting in and out
- // Supports mobility
- // Safely return to bed
- // Security



Individually adjustable comfort positions support residents' sensory training and allow them to participate in daily activities.

Ergonomic and functional lying surface for variable setting of care and comfort positions





Lying independently – safely and as comfortably as possible

The floorline® 9.5|80 offers an excellent range of features and is not only suitable for use as a lie-low bed, but also as a flexible and convenient nursing care bed for the care givers' specific needs. Residents with restless motor functions benefit in particular from the safety provided by the floorline® 9.5|80: in the low sleeping and resting position of only 9.5 cm above the floor, the risk of injury due to accidentally sliding out of the bed is nearly eliminated.

In addition, the barrier-free ultra-low position with lowered side rails makes it possible for many residents, should they happen to slide out, to get back into bed on their own.

Another feature that offers residents and care givers the best possible comfort and safety is the bed's ability to be lowered from 80 down to 9.5 cm, with a preset stop at the entry height of 37.5 cm as well as safety stops at 27.5 and 15.5 cm.

The floorline® 9.5|80 is available in various designs and with different end panels and also provides clearance for lifters and other aids to help residents out of bed.

The integrated telescoping bed and side rail extension to 220 cm is pulled out in no time.



Type 31



Type 32



Technical advantages at a glance:

- // Height adjustment: 9.5–80 cm
- // High lateral stability
- // Safe capacity: 220 kg
- // Outer dimensions: 102 × 240 cm
- // Smooth-rolling double-swivel castors (ø 100 mm) with integrated wall protection system
- // New, ergonomic lying surface
- // Safety via stops at 37.5, 27.5 and 15.5 cm
- // Integrated bed extension including side rail extension

.bock[®] ///



Would you like more information?

Additional information on the floorline[®] 9.5|80 can easily be found on the internet. Just scan the QR code shown and you will be automatically directed to our website.

Would you rather seek the advice of a professional? Then you can also reach us by telephone or email. Our enthusiastic team will be glad to help you!

Hermann Bock GmbH
Nickelstraße 12
33415 Verl
Germany

Phone: +49 5246 9205-0
Fax: +49 5246 9205-25
e-Mail: info@bock.net
Internet: www.bock.net